

## OPLC Ordering Information

Part #	Description
<b>Personal OPLC Instruments</b>	
<b>POSU 50</b>	
200001	<i>Purification unit supplied with R20 FEW rinsing cassette, A5, and A20 1 Channel FEW cassettes including SK003 starter kit.</i>
	Compatible with any HPLC system. 110/240 V, 50-60, 160VA
<b>Development Cassettes and PTFE sheets</b>	
200902	A-5 K7 cassette for 5x20 cm aluminium backed OPLC layers 200 µm thick
100900	A-20 K7 cassette for 20x20 cm aluminium backed OPLC layers 200 µm thick
100901	R-20 K7 cassette for rinsing and bidirectional separation on 20x20 cm aluminum backed OPLC layers 200 µm thick
100904S	G-20 SK7 cassette for 20x20 cm glass backed analytical OPLC layers 200 µm thick, with s.s. base plate
100906S	G-20 SK7 cassette for 20x20cm glass backed semi-preparative OPLC layers 500 µm thick, with s.s. base plate
100909S	G-20 SK7 cassette for 20x20cm glass backed preparative OPLC layers 1000 µm thick, with s.s. base plate
100908	PTFE sheets (upper and bottom) for R-20 rinsing K7 cassette
100910	PTFE sheets (upper and bottom) for A-20 K7 and G-20 SK7 type cassettes
100912	PTFE sheets (upper and bottom) for A-5 K7 cassette
<b>Development Cassettes and PTFE sheets, 1 Ch FEW</b>	
200900	A-20 FEW K7 for 20x20 cm aluminium backed OPLC layers 200 µm thick
200901	R-20 FEW K7 for rinsing and bidirectional separation on 20x20 cm aluminium backed OPLC layers 200 µm thick
200904S	G-20 FEW SK7 cassette for 20x20 cm glass backed analytical OPLC layers 200 µm thick, with s.s. base plate
200906S	G-20 FEW SK7 cassette for 20x20cm glass backed semi-preparative OPLC layers 500 µm thick, with s.s. base plate
200909S	G-20 FEW SK7 cassette for 20x20cm glass backed preparative OPLC layers 1000 µm thick, with s.s. base plate
200908	PTFE sheets (upper and bottom) for R-20 FEW rinsing K7 cassette
200910	PTFE sheets (upper and bottom) for A-20 FEW K7 and G-20 FEW SK7 type cassettes
<b>Starter kits</b>	
SK 001	Starter kit for A-20 type 20x20 cm K7 cassette, including:
	5 pieces, 20x20 cm Aluminium backed, 5 µm SiO <sub>2</sub> , F254, 200 µm thick LA 001 OPLC Layers
	5 pieces, 20x20 cm Aluminium backed, 11 µm SiO <sub>2</sub> , F254, 200 µm thick LA 011 OPLC Layers
SK 003	Starter kit for A-5 type 5x20 cm K7 cassette including:
	10 pieces 5x20 cm Aluminium backed, 5 µm SiO <sub>2</sub> , F254, 200 µm thick LA 003 OPLC Layers
	5 pieces, 5x20 cm, Aluminium backed, 11 µm SiO <sub>2</sub> , F254, 200 µm thick LA 013 OPLC Layers
LA 001	<b>OPLC Layers</b>
	<i>Aluminium backed, 5 µm particle size, 200 µm thick OPLC Layers</i>
	Pack of 20, 20x20cm, SiO <sub>2</sub> F254 (for A-20 Cassette) for analytical separations
LA 1001	Pack of 20, 20x20cm, SiO <sub>2</sub> (for A-20 Cassette) for analytical separations
LA 003	Pack of 40, 5x20cm, SiO <sub>2</sub> F254 (for A-5 Cassette) for analytical separations
LA 1003	Pack of 40, 5x20cm, SiO <sub>2</sub> (for A-5 Cassette) for analytical separations
LA 011	<i>Aluminium backed, 11 µm particle size, 200 µm thick OPLC Layers</i>
	Pack of 20, 20x20cm SiO <sub>2</sub> F254 (for A-20 Cassette) for analytical separations
	Pack of 20, 20x20cm SiO <sub>2</sub> (for A-20 Cassette) for analytical separations
	Pack of 40, 5x20cm SiO <sub>2</sub> F254 (for A-5 Cassette) for analytical separations
	Pack of 40, 5x20cm SiO <sub>2</sub> (for A-5 Cassette) for analytical separations
	Pack of 15, 20x20cm RP18 F254 (for A-20 Cassette) for analytical separations
	Pack of 25, 5x20cm RP18 F254 (for A-5 Cassette) for analytical separations
	Pack of 15, 20x20cm NH <sub>2</sub> F254 (for A-20 Cassette) for analytical separations
	Pack of 20 20x20 cm NH <sub>2</sub> F254 (for A-5 cassette) for analytical separations

## OPLC Ordering Information

Part #	Description
LA C11	Pack of 20, 20x20cm Al <sub>2</sub> O <sub>3</sub> 60 neutral F254 (for A-20 Cassette) for analytical separations
LA C13	Pack of 40, 5x20cm Al <sub>2</sub> O <sub>3</sub> 60 neutral F254 (for A-5 Cassette) for analytical separations
	<i>Glass backed, 11 μm particle size, 200 μm thick OPLC Layers</i>
LG 011	Pack of 20, 20x20cm SiO <sub>2</sub> F254 (for G-20 0.2mm Cassette) for analytical separations
LG 011-02	Pack of 10, 20x20cm SiO <sub>2</sub> F254 (for G-20 0.2mm Cassette) for analytical separations
	<i>Glass backed, 25 μm particle size, 500 μm thick OPLC Layers</i>
LPG 001	Pack of 15, 20x20 cm, SiO <sub>2</sub> F254 (for G-20 0.5mm Cassette) for semi-preparative separations
LPG 001-05	Pack of 10, 20x20 cm, SiO <sub>2</sub> F254 (for G-20 0.5mm Cassette) for semi-preparative separations
LPG 1001	Pack of 15, 20x20 cm, SiO <sub>2</sub> (for G-20 0.5mm Cassette) for semi-preparative separations
LPG 1001-05	Pack of 10, 20x20 cm, SiO <sub>2</sub> (for G-20 0.5mm Cassette) for semi-preparative separations